

IMAGING PROTOCOL

for digital planning in orthognathic surgery
CT / CBCT acquisition and dental data only

ACCURACY · REPRODUCIBILITY · PREDICTABILITY




The quality of digital planning depends directly on the quality of the imaging acquisition and the dental data. Following the protocol = surgical predictability.

1 CT / CBCT ACQUISITION

GENERAL PRINCIPLES

- Scan in an upright position (seated / standing)
- Natural Head Position (NHP)
- Mandible in centric relation (CR)
- No chin support
- No forehead support / bands across the forehead
- Patient in slight disocclusion during the scan



RECOMMENDED


Stabilise the head with posterior support and a light forehead band only if necessary.

DO NOT USE

Chin support (it deforms the soft tissues) or any pressure on the maxilla.

IMPORTANT — METAL ARTEFACT REDUCTION


Enable the metal-artefact-reduction protocol during CT / CBCT acquisition — it reduces streaks and shadows from metallic materials (braces, fillings, implants).



2 RECOMMENDED CT / CBCT PARAMETERS

PARAMETER	CT (SPIRAL)	CBCT (CONE-BEAM)
FOV (field of view)	Full head (min. 16–22 cm)	Full head (min. 16–22 cm)
Matrix	≥ 512 × 512	≥ 512 × 512
Slice thickness (Z)	0.5–1.0 mm (reconstruction)	0.3–0.6 mm
In-plane reconstruction (X,Y)	0.5–0.7 mm	0.3–0.4 mm
Orientation	Axial, preferably	Axial, preferably
Contrast agent	NO	NO
Scan time	Per device protocol	20–40 sec
Table feed (CT spiral)	< 3 mm / rotation	—
Metal-artefact-reduction protocol	ENABLED	ENABLED

2.1 — SCAN VERIFICATION (SCOUT)



Wax wafer: a wax occlusal record of the CR position — use it to keep the mandible in centric relation during the scan, if necessary.

2.2 — OCCLUSION IN CR (CENTRIC RELATION)

✓ CORRECT — IN CR

- Uniform bilateral contacts
- Relaxed mandible
- Condyles centric



✗ INCORRECT — OTHER POSITION

- Anterior / lateral shift
- Non-uniform contacts
- Condyles not centric



3 DENTAL DATA (MANDATORY)

CBCT alone does not provide precise occlusal information on its own — dental data completes the picture.

! THREE SCANS ARE REQUIRED

We always need a dedicated scan of each arch plus the bite — the CT / CBCT alone is not enough:

- **Upper jaw** — full maxillary arch
- **Lower jaw** — full mandibular arch
- **Occlusion** — both arches together, in centric relation (CR)

3.1 — OPTIONS FOR CAPTURING DENTAL DATA



1. Intraoral scan (preferred)

- High accuracy · STL directly usable · reproducible



2. Plaster models

- Scanned with an optical scanner or lab CBCT · includes the CR occlusion

3.2 — FILE FORMATS: DENTAL DATA

Intraoral scan	.STL (binary / ASCII) · .PLY · .OBJ
Scanned models (optical / CBCT)	.STL · .PLY · .OBJ
Occlusion in CR	.STL — both arches, same relation

3.3 — FILE FORMATS: CT / CBCT

Imaging data	.dcm · .DICOM (uncompressed)
Reconstructions / surface	.dcm · .DICOM
Volume export (optional)	DICOMDIR (.dcm series)

GENERAL RECOMMENDATIONS

- Always provide the complete data set (all DICOM files)
- Avoid lossy compression on export
- Include patient and orientation information
- Communicate the CR position and the type of record used

FOR OPTIMAL RESULTS

- Always follow the positioning and scan parameters
- Provide complete, correct dental data
- Share the relevant details with the planning team

SIMPLIFIED WORKFLOW



KEY MESSAGE

Planning quality depends on two essentials:
1. Correct CT / CBCT acquisition per protocol
2. Complete, precise dental data in CR