

IMAGING PROTOCOL

for patient-specific implants

CT / CBCT acquisition — the defect & the healthy mirror side

ACCURACY · SYMMETRY · RECONSTRUCTION




A patient-specific implant restores a skeletal defect. Its fit depends directly on the CT — the defect margins and the healthy mirror side must be captured cleanly.

1 CT / CBCT ACQUISITION

GENERAL PRINCIPLES

- Capture the full defect region (cranial · orbital · midface · mandibular)
- Include the contralateral healthy side — for mirroring
- Natural Head Position — stable and reproducible
- Patient still — avoid any motion during the scan
- Metal-artefact-reduction protocol enabled



Stable head position — defect & mirror side in the field of view

RECOMMENDED

Capture the defect and the intact opposite side in the same scan, so the implant can be mirrored from the healthy anatomy.

IMPORTANT — METAL ARTEFACT REDUCTION


Enable the metal-artefact-reduction protocol — existing hardware, plates or dental work near the defect otherwise cast streaks over the reconstruction surface.



2 RECOMMENDED CT / CBCT PARAMETERS

PARAMETER	CT (SPIRAL)	CBCT (CONE-BEAM)
FOV (field of view)	Defect + contralateral side	Defect + contralateral side
Matrix	≥ 512 × 512	≥ 512 × 512
Slice thickness (Z)	0.5–1.0 mm (reconstruction)	0.3–0.6 mm
In-plane reconstruction (X,Y)	0.5–0.7 mm	0.3–0.4 mm
Orientation	Axial, preferably	Axial, preferably
Contrast agent	NO	NO
Scan time	Per device protocol	20–40 sec
Table feed (CT spiral)	< 3 mm / rotation	—
Metal-artefact-reduction protocol	ENABLED	ENABLED

2.1 — SCAN VERIFICATION (SCOUT)




Defect covered Mirror side included No motion Stable position

Tell us the defect site — cranial, orbital, midface or mandibular. Capture the defect and the mirror side at full resolution so the implant can be designed and mirrored accurately.

3 DATA YOU PROVIDE

NO DENTAL DATA REQUIRED

For cranial, orbital and midface implants, the design works from the bony defect. **Intraoral scans and occlusal records are not needed.**



The implant is designed to the defect and mirrored from the healthy side

3.1 — WHAT TO PROVIDE

PROVIDE

- One CT or CBCT covering the defect **and the mirror side**
- Metal-artefact-reduction enabled
- The complete, uncompressed DICOM data set
- The defect site and the goal (restore / augment / support)

MANDIBULAR / TMJ IMPLANTS

Only for implants that restore the mandible or the TMJ: scan with the condyles seated in CR, and add a dental scan if the bite is involved.

3.2 — FILE FORMATS: CT / CBCT

Imaging data	.dcm · .DICOM (uncompressed)
Reconstructions / surface	.dcm · .DICOM
Volume export (optional)	DICOMDIR (.dcm series)


GENERAL RECOMMENDATIONS

- Always provide the complete data set (all DICOM files)
- Avoid lossy compression on export
- Include patient and orientation information
- Tell us the defect site and the reconstruction goal

FOR OPTIMAL RESULTS

- Follow the positioning and scan parameters
- Capture the defect and the healthy mirror side
- Share the defect site and the reconstruction goal

SIMPLIFIED WORKFLOW



KEY MESSAGE

Implant accuracy depends on one essential:

A clean CT / CBCT of the defect and the healthy mirror side. Dental data only for mandibular / TMJ implants.